

## Typical OR Equipment

- Upon entering any OR your first priority should be to make sure the room has all the equipment required for procurement - and arranging said equipment appropriately.
- **\*\* If hospital has a slush machine you are utilizing START IT NOW \*\***
- In addition to the supplies below you will need tables and IV poles - make sure there are enough based on the organs being procured

### OR Table/Bed -

- Make sure your OR bed has all parts (head and foot are removable) and is orientated appropriately. A sheet or back-table cover should be used to cover cushions on bed.
- Make sure bed is locked and has power - clip remote to head of bed for easy access.
- Arrange EKG leads (+ stickers) on head of bed with left leads hanging off left side and right leads hanging off right side. Place folded blanket on top to hold them in place and create a 'pillow' for donor.
- Bovie grounding pads can also be gathered placed to the side.
- Make a shoulder roll to place under donor (check donor weight to judge size of roll required)
- Pull additional 2-3 blankets at this time - they will come in handy for covering patients' legs and soaking up blood/fluids on floor.
- Pull all prep supplies for DCD - if BD case consult scrub and circulator to see what supplies they may already have open/pulled and their expectations on who will prep.



### Bovie/Cautery/Bipolar

- For BD abdominal case: two bovie are required. If abdominal and thoracic: 3 bovie
- Generally, our surgeons like to have the bovie set on 80/80 at start of cases.
- Pancreas: a Liga-Sure will most likely be used - ask circulator to have it ready in room
- For DCD: no bovie required (\*If heart team is using ECMO method and restarting heart, plan as if case is a BD donor)



### Suction/Neptune

- One completely empty suction is required for our cases. If 1st Neptune is not completely empty or if large number of organs are being procured have a second Neptune close-by.
- Connect manifold (A) in lower chamber position - connect all four suction lines to bottom chamber.



## Bronch

- Required for all lung cases.



## Crash Cart with internal paddles

- Required for all thoracic cases - don't open just have in room or right outside door.



'Standard' OR



### Scrub Supplies

- Once your room is set up - check on your scrub to make sure they have all the necessary instrumentation
- Key instruments:
  - 'LOOP' tray or Major set (many hospitals have a "LOOP set" but the contents vary to make sure they have all the essentials)
    - Multiple Metz
    - Multiple Mayo scissors (curved and straight)
    - Multiple Tonsils
    - Multiple right angles
    - Multiple Hemostats (aka Snaps)
    - Multiple Debakeys (aka pickups, aka forceps)
    - Needle drivers (fine tip if possible)
    - 2-3 Knife handles
    - Kelley's
    - Kocher's
  - 4 penetrating towel clamps
  - Vascular Clamp set (if the LOOP set does not have cross-clamps)
  - Sternal Retractor

- Sternal Saw and blade
  - Vascular dilators
- Thoracic Case:
  - Open heart set useful.
    - Ask teams if they need anything specific
- Other supplies:
  - 4 suction tubing
  - 2 (or 3) bovie
  - Asepto Syringe
  - Multiple basins
  - 0-0 and 2-0 silk ties
  - Bone wax (BD only)
  - Umbilical tape
  - Kidney biopsies: 5-0 prolene suture x2
  - Heart: 5-0 silk pop-offs multi-pack